

User Manual

COBRA



A.R.C. Laser GmbH

CE 0123

Bessemerstr. 14
90411 Nürnberg

Tel.: +49 911 21779-0
Fax: +49 911 21779-99

info@arclaser.de
www.arclaser.de

Please contact for service:



A.R.C. Laser GmbH Headquarters:

A.R.C. Laser GmbH
Bessemerstr. 14
D-90411 Nürnberg

Tel.: +49 911 21779-0
Fax: +49 911 21779-99
service@arclaser.de

Copyright © A.R.C. Laser GmbH. All rights reserved.

This document contains confidential and copyrighted information and may only be used with the devices described herein. Reproduction is prohibited without written permission from A.R.C. Laser GmbH.

1	INTRODUCTION.....	5
1.1	Marking and Symbols.....	5
1.2	Intended purpose.....	6
1.3	Characterization of the User	7
1.4	Characterization of the Patients	8
1.5	Benefit Claims	9
1.6	Theory and Technical Set-up.....	9
2	Transport and Storage	10
2.1	Shipping and unpacking the device.....	10
2.2	Return Shipment	11
3	Set up and installation	12
3.1	Installation Site	12
3.2	Room requirement.....	12
3.3	Electrical Connection	13
4	Safety Information	14
4.1	General.....	14
4.2	Eye Safety.....	15
4.3	Electrical Protection	15
4.4	Explosion and Fire Hazard.....	16
4.5	Protection against Undesired Radiation	16
4.6	NOHD Safety Distance.....	16
4.7	CE-Regulations	17
4.8	RoHS3- Regulation	17
4.9	Protective Housing.....	17
4.10	Connectors and switches	18
4.11	External Interlock Switch	19
4.12	Safety-Shutter/Aiming Beam	19
4.13	Manual Reset	19
4.14	Reset by Power Failure	19
4.15	Labels and Markings	20
4.16	Operating Conditions:.....	22
4.17	Electromagnetic Compatibility.....	22
5	Advice foR users.....	23
5.1	Technical Instruction.....	23
5.2	Laser Safety Training	23
5.3	Medical Instruction	23
5.4	Medical Device Book.....	23
5.5	Device Parts and Accessories	24
6	Operation.....	27
6.1	Start the Cobra.....	28
6.2	Operation of CITO 532	32
6.3	Operation of the Q-Las.....	36
6.4	Switch Off your Cobra	39
6.5	Laser-Stop	39
7	SPECIFICATIONS	40
7.1	General.....	40
7.2	Laser data.....	40
7.3	Electrical Connection Data	41
7.4	Classification	41
8	Service	42

8.1	Introduction	42
8.2	Safety Check (STK).....	42
8.3	Care by the user	43
8.4	Slit lamp	43
8.5	Troubleshooting.....	45
8.6	Disposal.....	49
9	CUSTOMER service.....	50
9.1	Warranty Information.....	50
9.2	Warranty, Shipment, Packing	50
9.3	Sales and Service Information.....	50
10	GUIDELINES And MANUFACTURERS DECLARATION	51
10.1	Electromagnetic Emissions	51
10.2	Electromagnetic Immunity (1)	52
10.3	Electromagnetic Immunity (2)	53
10.4	Recommended separation distances between portable and mobile RF telecommunications equipment and the laser	54

1 INTRODUCTION

We thank you for having chosen our Cobra laser system and wish you every success in using this versatile laser system.

Cobra Laser system generates a laser beam of high intensity, which can cause injuries if handled improperly. The user manual should therefore be read carefully before starting up the device. If you have any further questions regarding the safety or the use of the device or regarding the laser and laser radiation, please contact A.R.C. Laser GmbH.

1.1 Marking and Symbols



The "**General warning sign**" (ISO 7010 W001) indicates general risk. It encourages the user to take care regarding the hazard specified by the supplementary sign. Before carrying out any further works at parts with such marks, please read the user manual or contact your local dealer or contact directly the A.R.C. Laser GmbH service department.



The sign "Warning; Laser beam" (ISO 7010 W004) indicates laser hazards. It encourages the user to take care to avoid exposure to a laser beam. The laser beam of this device **may not be visible** to the naked eye nor through protecting goggles.

1.2 Intended purpose

The device COBRA combines Q LAS (Nd:YAG Laser) and CITO 532 (SLT Laser) at one slit lamp. The major function of COBRA is switching back and forth between the two lasers CITO 532 and Q LAS. Therefore, no further purpose is intended in addition to CITO 532 and Q LAS.

1.2.1 General Purpose

Irradiation of trabecular meshwork by laser light, laser light is acting on patient's tissue without destroying tissue [CITO 532]. Tissue disruption by plasma shockwave in the human eye [Q LAS].

1.2.2 Medical Purpose

The COBRA itself has only an indirect medical purpose. The purpose of the COBRA is switching back and forth between the two lasers CITO 532 and Q-LAS. The general purpose in chapter 2.1 depends on the used laser.

In the following chapter the medical application, indications, contraindications and side and adverse effects for the devices of CITO 532 and Q-LAS are shown.

1.2.3 Indications and contraindications

Indirect, depends on the used laser CITO 532 or Q-LAS.

Device	Medical Application indicated	Indication from clinical data	Side and adverse effects	Contraindications
CITO 532	Open angle glaucoma/ Selective laser trabeculoplasty	-Ocular hypertension in combination with an open chamber angle (POAG or OAG)	-discomfort (pain) -IOP rise (non-permanent) -vision decrease (non-permanent) -inflammation, keratitis -photophobia -hyperemia -development of synechia -corneal oedema, haze, burns -uveitis -retinal detachment -cataract -blurred vision (not permanent) -headache (not permanent)	-Closed chamber angle -congenital glaucoma -inflammatory/ uveitic glaucoma -poor visualization of trabecular meshwork
Q-LAS	Iridotomy	prevention or treatment of high intraocular pressure related to angle closure glaucoma	-IOP rise (non-permanent) -inflammation -bleeding -change of endothelial cell density -cystoid macular edema -retinal complications -discomfort (mild/moderate)	-open angle glaucoma -congenital glaucoma -inflammatory/ uveitis glaucoma -acute intraocular inflammation -acute reduced optical transparency of the cornea

			-pain -scotoma	
Q-LAS	Capsulotomy	reduced visual acuity related to PCO	-inflammation -change of endothelial cell density -cystoid macular edema -retinal detachment -IOP rise -bleeding -lens damage / pitting lens	-acute intraocular inflammation -acute reduced optical transparency of the cornea

ATTENTION

The device may only be operated by specially trained personnel who are experts in the medical effects and possible dangers of the device. You must have the necessary skills to use the laser in accordance with this instruction manual.

When not in use, the device should be protected against unauthorized use.

1.3 Characterization of the User

Laser systems may only be operated by doctors and healthcare professionals. Each user must be trained in the operation of the laser by a specialist authorized by A.R.C. Laser or by someone trained by an authorized specialist.

1.3.1 Ophthalmologist

Typical job title	Ophthalmologist
Provided education	Medical degree, Ophthalmology specialist
Demographic particularity	Not necessarily native speakers, speaks language of the organization at least at B2 level
Provided work experience	Training with experienced doctor is recommended
Typical work environment	Treatment room (s. Chapter 4.4)
Typical work	Treatment of the patient
Provided training	No training provided. Training by a A.R.C. Laser specialist or person trained by a A.R.C. Laser specialist is recommended.

1.3.2 Medical assistant

Typical job title	Medical assistant
Provided education	Vocational training
Demographic particularity	Not necessarily native speakers, speaks language of the organization at least at B2 level
Provided work experience	N/A
Typical work environment	Treatment room (s. Chapter 4.4)
Typical work	Cleaning the device, disinfect the device, prepare and assist the surgery
Provided training	No training provided. Training by a A.R.C. Laser specialist or person trained by a A.R.C. Laser specialist is recommended.

1.3.3 A.R.C. Laser Service Employee

Typical job title	A.R.C. Laser Service Employee, medical technician
Provided education	Vocational training, Training by experienced A.R.C. employees for COBRA
Demographic particularity	Not necessarily native speakers, speaks language of the organization at least at B2 level
Provided work experience	N/A
Typical work environment	Treatment room (s. Chapter 4.4) service department of A.R.C.
Typical work	Installation, Servicing, Safety Check
Provided training	Training for COBRA

1.4 Characterization of the Patients

The *COBRA* laser device is used in the field of Ophthalmology. In general, the *COBRA* can be used for all patients, where the indications of the *CITO 532* or *Q-LAS* meet. Please note the contraindications for the use of both devices.

Exceptions may be when the patient is not able to sit in front of the slit lamp, open the eye and keep it in the same position for the duration of the treatment.

Patients under 18 years are excluded from the treatments.

The prevalence of the diseases predominantly increases with age.

treatment-specific requirements:

SLT: The treatment can be used for the indications (P)OAG and OHT. The trabecular meshwork must be visible for treatment.

SLT is also a good treatment option for pregnant and breastfeeding women, as medication is contraindicated in this population.

Iridotomy: The treatment can be used for the indication ACG or as preventive action to avoid ACG.

Capsulotomy: The treatment can be used for the indication PCO.

1.5 Benefit Claims

A combination laser device is more cost-effective as only one table and one slit lamp is needed. In addition, a combination device is more space-saving than two individual laser devices and offers the physician a more diverse range of treatment-options than a single laser.

1.6 Theory and Technical Set-up

Cobra is a dual platform laser device consisting of the CITO 532 and Q-Las, featuring full performance capabilities of both devices. The laser beam path is concurring with the observation path of the slit lamp. Using the aiming beam, the desired location of the treatment can be precisely positioned.

1.6.1 Theory and Technical Set-up CITO 532

The CITO 532 is a diode-pumped laser with a wavelength of 532 nm. The radiation in the laser head is generated with the help of a pump diode. The laser beam that leaves the device has a small diameter; this is 400 µm at the point action in the patient's eye.

The wavelength of the emitted laser light depends on the laser medium. In our case, infrared laser radiation is generated in the device. The wavelength of this infrared laser radiation is now halved with the help of a special **crystal** and generates a wavelength of 532 nm. The laser beam is converted into short, green laser pulses with high energy, which leaves the device through the laser beam path.

With the CITO 532, the laser beam path is overlaid with the observation beam path of the slit lamp. With the help of a target beam, the desired treatment position can be precisely determined in advance.

1.6.2 Theory and Technical Set-up Q-Las

The Q-Las is a Nd: YAG laser with a wavelength of 1064 nm. The radiation in the laser head is generated with the help of a flash lamp or a diode (depending on the purchased Cobra variant). The laser beam that leaves the device has a small diameter; this is less than 10 µm at the place of its effect in the patient's eye.


The wavelength of the emitted laser light depends on the laser medium. In our case, infrared laser radiation is generated. Attenuating optics enable the energy that should leave the device to be set. Internal control mechanisms compares the actual with the target state for the selected amount of energy.

The focus of the laser can be shifted in 2 steps using a dial on the laser unit.

2 TRANSPORT AND STORAGE

We at A.R.C. Laser GmbH will make sure that the device is packed and transported with the greatest possible care.

Before unpacking the laser system, please check the packaging for damage and report any damage immediately to the shipping agent and A.R.C. Laser GmbH. Only remove the packaging in the presence of a representative of the carrier. Make a list of the damaged parts and have this list signed by the courier.



The device must be transported at temperatures between 0 °C and +50 °C. The air pressure during transport must be between 1080 hPa and 750 hPa. During storage, a temperature range of 5 °C to 40 °C must be maintained. The environment/air must be dry and clean. The relative humidity during transport and storage must be between 5 % and 85 %.

ATTENTION

If the laser is transported or stored at a temperature below 5 °C, it can be damaged when starting. Unpack the laser and leave it for at least half a day (12 hours) at normal room temperature so that the system reaches room temperature.


2.1 Shipping and unpacking the device

Unpacking and installation of the device must be carried out by an authorized technician or authorized representative of A.R.C. Laser GmbH. After the installation and the correct connection of the device, the technician will put the device into operation and explain the function of the device. All functions and safety procedures are discussed and demonstrated here.


2.2 Return Shipment

The device will be shipped in a specially developed transport packaging. Please keep this. If you ship the device back to A.R.C. Laser GmbH, be sure to use the intended transport packaging.

The transport packaging is designed so that the device fits snugly into the foam parts.



All accessories are stowed in an extra padding inside the transport packaging. If all foam parts are in the right place and the box is closed, the device is ready for shipping.



3 SET UP AND INSTALLATION

3.1 Installation Site

Before the device is delivered, it must be ensured that the laser can be set up in a suitable location.

The Cobra Laser should be operated in an easily accessible place. The laser should not be operated near a heater, since air cooling works best when the ambient temperature is not higher than 28 °C. Higher temperatures can cause the device to switch off due to overheating. A place of installation in direct sunlight can cause an early shutdown and must be avoided in any case.

If the temperature is too low (below 18 °C), the device cannot be started to prevent condensation on the internal optics; this could result in permanent damage to the laser.

- The device should be set up so that the laser beam is not aimed directly at a door, window or reflective material.
- The device should be operated in a dust-free room. There should be no carpets on the floor or the walls.
- When the device is not in use, it should be covered to prevent contamination. For this purpose, a suitable dust cover is included.

The wiring must be installed so that there are no tripping hazards or other hazards

All control elements must also be free and easily accessible. The air humidity is monitored internally and must be below 75 %.

3.2 Room requirement

The legislation imposes the following requirements to any room in which a class 3B laser according to IEC 60825-1 is operated.

3.2.1 Marking Access Points

All entrances must be clearly marked so that unintentional or unauthorized entry, which can lead to a hazard, is avoided.

- Laser warning signs (triangle with laser symbol) and the wavelength labeling must be attached to each access door.
- A warning lamp must be attached above each access door. This must always light up when the laser is in operation.
- Unintentional entry into the room without safety glasses is avoided.
- The laser safety goggles must be easily accessible at the entrance.

3.2.2 Window Shielding

It must be ensured that no laser radiation can escape the room. In particular, the windows must be secured with suitable materials. If you have any questions, your contact person at A.R.C. Laser GmbH is available to assist you.

3.2.3 Reflecting Surfaces

To avoid danger from reflected direct and scattered radiation, there must be no reflective surfaces in the room.

These can be:

- Mirrors
- Pictures behind glass
- Chrome surfaces
- Windows

These surfaces must be removed or suspended or matted. In the area of laser use, only matted, non-reflective and non-flammable instruments and materials should be used.

3.3 Electrical Connection

The laser may only be operated on an earthed wall socket and can be operated from 100 V - 240 V (50 Hz/60 Hz) AC voltage.

The power connection is made via the table column and the supplied IEC connector.

ATTENTION

Make sure that the plug is accessible at any time so that the laser can be disconnected from the mains after use.

4 SAFETY INFORMATION

4.1 General

The Cobra Laser system is a precision instrument for medical applications. The system has been carefully developed and tested through intensive testing before shipping. In order to offer you and the operating personnel all possible protection, we recommend that you read this section of the operating instructions carefully.

The Cobra Laser is a class 3B laser according to EN 60601-2-22 or rather EN 60825-1.

Class 3B in the standard describes high-energy lasers and therefore special measures must be taken before commissioning in order to ensure safe and trouble-free working with the device. It is particularly important to protect the eyes and skin of the operator, the patient and third parties. Laser safety glasses must be used for eye protection.

The following explanations are not exhaustive. All users of laser devices should enclose applicable legal regulations and provisions with the device and inform the staff accordingly.

If the device is operated outside Germany, provisions of the American National Standard Office ANSI Z136.3-2018 "American National Standard for the Safe Use of Lasers in Health Care Facilities" and ANSI Z136.1-2014 "American National Standard for the Use of Lasers" must be respected.

This manual is limited to the operation, maintenance, and control of the device. The manual is not a guide for the treatment of diseases that can be remedied by laser.

With regards to the devices supplied, such as slit lamps or instrument tables, the safety, operating and maintenance instructions in the corresponding manuals must be observed. With regards to the instrument table, reference is only made here to the danger from improper use as a seat or storage area. When operating the height adjustment, make sure that no one can be harmed.

A.R.C. Laser GmbH cannot be held responsible for damage or damage resulting from improper use.

The warranty for the device expires if the laser has been opened (even partially), modified or repaired by unqualified personnel.

4.2 Eye Safety

As a safety measure against direct or indirect laser radiation, it is necessary that everyone in the room wear laser safety glasses. The treating doctor is protected from radiation by the laser protection integrated (eye safety filter) into the slit lamp. Eye protection appropriate to the laser must be guaranteed for the patient.

When using the laser, only laser safety glasses that are designed for the respective wavelength and on which the CE mark and the appropriate protection class are noted, may be used.

The protection level (LB) for laser safety glasses is determined by the standard EN 207 (personal eye protection). The optical density (OD) describes the attenuation of the radiation through the glasses.

There are different types of laser safety glasses that also allow spectacle wearers to wear a tight goggle on top of their own which are protective to all sides. In particular, the scattered radiation that does not directly come from the front of the eye may present risks resulting from internal reflections on the glasses. Therefore, we encourage you to wear laser safety glasses which also guarantee a full protection to the side parts.

When using the CITO 532, you must use safety goggles for **532 nm** and provide safety class **DIR LB5 (OD 5+)**.

By using the Q-Las, you must use safety goggles for **1064 nm** and provide safety class **DIR LB6 (OD6+)**.

Alternatively, there are also laser protective goggles that are designed for both wavelengths.

A.R.C. Laser GmbH offers laser safety glasses for the use of the Cobra. See chapter 5.5.3 for detailed information.

ATTENTION

Never look directly at the laser beam or the light reflected by the laser beam.

This will cause serious damage to your retina.

4.3 Electrical Protection

Never remove the housing parts of the laser. Work on the device may only be carried out by authorized A.R.C. Laser GmbH service personnel. If the mains have to be replaced, it must be carried out by authorized service personnel only.

The room in which the laser is operated should be kept dry. In the event that cleaning is necessary, make sure that the floor is dry before starting up the device.

ATTENTION

Never use the laser if you notice any visible defect on the device.

Never use the device if you notice any visible defect on the power plug, wires or supply lines are exposed.

4.4 Explosion and Fire Hazard

Never work with the laser near flammable anesthetics, flammable solutions or other flammable materials. Flammable plastic parts or paper parts in particular should also be removed from the working area of the laser. There is a fire or explosion hazard if the laser is used in the presence of flammable materials, solutions, or gases or an oxygen-enriched environment.

4.5 Protection against Undesired Radiation

If a footswitch is used, it should never be outside the area of the attending doctor. It is forbidden for anyone other than the attending doctor to trigger the footswitch.

Especially in operating theaters where there are multiple footswitches, it is important to ensure that the laser footswitch is close to the laser device.

4.6 NOHD Safety Distance

The NOHD (Nominal Ocular Hazard Distance) is the distance at which the irradiance is equal to the exposure limit value of the cornea of the eye. The NOHD identifies the danger area within which there is a risk of damage to the health of the eye if you look directly into the laser beam without protection.

The NOHD is calculated according to EN 60825-1 and taking into account the permitted power fluctuations (+/- 20%) according to EN 60601-2-22 using the following equation:

$$NOHD = \frac{\sqrt{\frac{4P}{MPR * \pi} - Diameter\ Beam\ Bundle}}{Beam\ Divergence}$$

Calculation of the maximum permissible irradiation (MPR) with an exposure time of 100s

For the **CITO 532** applies:

Wavelength λ:	532 nm
NOHD:	3.7 m
Beam divergence $\alpha_{(full\ angle)}$:	11.7°
Permitted Maximum Radiation:	$5.34 \times 10^{-3} \text{ Jm}^{-2}$

For the **Q-Las** with flash lamp applies:

Wavelength λ:	1064 nm
NOHD:	3.6 m
Beam divergence $\alpha_{(full\ angle)}$:	16°
Permitted Maximum Radiation:	$2 \times 10^{-2} \text{ Jm}^{-2}$

For the Q-Las with diode the following NOHD applies:

Wavelength λ:	1064 nm
NOHD:	3.6 m
Beam divergence $\alpha_{(full\ angle)}$:	16°
Permitted Maximum Radiation:	$2 \times 10^{-2} \text{ Jm}^{-2}$

4.7 CE-Regulations

The laser system Cobra has been approved by the notified body in accordance with the European Regulation 2017/745/EU for medical devices. Accordingly, the device bears the CE mark **CE 0123**.

The device has been checked for electrical and mechanical safety. All parts we use comply with the CE regulations or have been tested for approval or suitability by the notified body.

Additional devices that you attach to the device require approval from an official test center. Changes to the device or interventions on your part will void the approval and warranty.

A device book and the test approval number are included with the device.

4.8 RoHS3- Regulation

Our company operates worldwide and regards the protection of the environment and natural resources as an entrepreneurial obligation. Based on individual tests, A.R.C. Laser GmbH confirms that, to the best of our knowledge, our products do not contain any substances in concentrations whose placing on the market is prohibited according to the applicable requirements of Directive 2015/863/EU (RoHS 3).

4.9 Protective Housing

The Cobra laser system has multiple protective housings. They prevent radiation of the lasers from escaping and protects users from touching live parts. This housing must not be removed.



1) Housing of the slit lamp PCL5,
the Q-LAS laser parts, and the
switching optics

2) Controller / Touchscreen

3) Housing of the CITO laser parts
and electrical components



4.10 Connectors and switches



The main switch under the table turns the entire system on and off.

To prevent a defect in the unit - caused by the power supply unit - the Cobra has two thermal fuses under the table. When the thermal fuses blow, they are pushed down and the red LED on the main switch goes out. As a result, no more voltage is passed on to the system and you cannot continue treatment. In this case please contact the service department of A.R.C. Laser GmbH or their local representative.

You can connect a warning light to the orange marked connector. ● This indicates to third parties whether the laser is currently in operation. The switch and wiring must be designed for at least 12 V and 100 mA.

The red marked connector is the interlock connector. ● You can replace it with a door interlock switch (see chapter 4.11).

There is a USB connection on the touch display. This is to be used for service purposes only! It is not possible to load or save files or programs via this port.

ATTENTION

To avoid the risk of electric shock, this device may only be connected to a supply network with a protective earth conductor.

4.11 External Interlock Switch

A door interlock switch is required by the accident prevention regulations. The device is equipped with an interlock connector as standard, which can be replaced by a door interlock switch. The laser switches off when the door is opened. In addition, the laser cannot be switched to the READY mode when the door is open. An error message appears on the control panel. When the door is closed, the error message disappears, and the laser can be switched READY again.



When installing a door interlock switch, note the following:

The switch and wiring must be provided for at least 12 V and 20 mA. The wires should end with a standard male connector. The choice of polarity is expedient in both variants. The socket for connecting the door switch can be found underneath the table. Insert the door switch connector there.



Make sure the socket is stuck firmly to prevent unexpected system interlock problems.

4.12 Safety-Shutter/Aiming Beam

The Cobra Laser has an internal safety shutter. This shutter is opened by pressing the READY button and then releases laser radiation. However, this can only be done if the laser has passed the internal tests and calibrations.

The aiming beam is only visible when the laser is in the READY mode. This aiming beam is a laser beam with very low power.

4.13 Manual Reset

If an error occurs, the system switches into STANDBY mode. You should then switch the laser off and on again. The restart should fix the error due to the automatic recalibration.

If the error still occurs, this can only be remedied by trained personnel. Please contact the service of A.R.C. Laser GmbH in this case.


4.14 Reset by Power Failure

If the device is disconnected from the power grid accidentally - e.g. due to a power failure - it usually restarts automatically. The system recalibrates and deletes all unsaved settings. If there is still an error, it can only be corrected by trained personnel. Please contact A.R.C. Laser GmbH service department should this occur.

4.15 Labels and Markings


The Cobra is provided with various warnings in accordance with European and worldwide guidelines. This is to prevent the laser user from being exposed to laser radiation due to carelessness.

4.15.1 Type Plate



- 1) Manufacturer
- 2) Laser name
- 3) Laser specifications
- 4) Electrical connection data
- 5) Serial number (SN)
- 6) Production date
- 7) GTIN (Global Trade Item Number)
- 8) UDI (Unique Device Identification = GTIN + SN + production date)
- 9) CE Mark
- 10) Medical Device
- 11) Follow instructions for use
- 12) Application part type BF
- 13) Do not dispose in household trash

4.15.2 Fuse Label



Only use the described fuses.

4.15.3 Modification Label




Shows the current device status based on the marked modification numbers.


4.15.4 Warning Label Laser



The "Laser Aperture" label marks the laser output. The other labels provide information about the type and intensity of the laser system.



4.15.5 Serial Number



The Cobra consists of several subsystems, each with its own serial number.

4.15.6 Laser-Stop



This symbol indicates the laser stop / emergency stop button (see chapter 6.5).

4.15.7 Follow instructions for use



On the device, you will again be reminded that the instructions for use must be observed.

4.16 Operating Conditions:

The medical laser Cobra is not suitable for use in connection with combustible gas mixtures of all kinds.

The device is not approved for operation at altitudes above 2,000 m above sea level. and only for an air pressure between 1080 hPa and 750 hPa.

The following environmental conditions must be met:

- Ambient temperature: 18 °C - 28 °C
- Relative humidity: < 75 %

4.17 Electromagnetic Compatibility

The system meets the EMC requirements according to EN 60601-1-2. Guidelines and the manufacturer's declaration are described in chapter 10.

ATTENTION

This device should not be operated immediately adjacent to or stacked with other devices.

This device should not be operated immediately adjacent to portable or mobile wireless communication devices or stacked with portable or mobile wireless communication devices.

For recommended separation distances from portable or mobile wireless communication devices, please refer to Chap. 10.3 (Electromagnetic Immunity - II) and 10.4 (Recommended separation distances).

RFID systems shall not be used at separation distances less than 0.15 m from the device. Otherwise, degradation of the performance of this device could result.

5 ADVICE FOR USERS

5.1 Technical Instruction

During the installation of the device, instruction is given by an A.R.C. Laser GmbH employee or an authorized representative.

This first instruction essentially relates to the technical use of the device. In addition, essential security-relevant points are dealt with. During the briefing, all persons working in the vicinity of the laser should be present.

After instruction, all instructed persons are noted in the device book, with one person being entered as the person responsible for safety. This person is later entitled to instruct other people on the device. These must also be noted in the device book.

5.2 Laser Safety Training

The Cobra laser is designed for medical users. It may only be used by personnel who have been instructed in its operation. The A.R.C. Laser GmbH recommends, in addition to the briefing, participation in seminars in which working with different laser systems is dealt with. In addition, instructions are given on laser safety and the use of lasers in general. It has also proven useful that people who do not work directly with the laser attend courses on laser safety.

Training for the accompanying staff is additionally accompanied by an instruction from an A.R.C. Laser GmbH employee or by an authorized representative when installing the device. During instruction, the use of laser safety glasses and laser safety will be specifically addressed.

A.R.C. Laser GmbH has a list of recommended courses as well as laser safety courses. These can be obtained from us at any time.

5.3 Medical Instruction

In the context of device instruction, only the general medical application is addressed. The A.R.C. Laser GmbH only gives recommendations for applications.

If necessary, there is the possibility to take part in a training course with an experienced doctor. Please contact your responsible administrator or A.R.C. Laser GmbH directly.

5.4 Medical Device Book


The medical device book is included with the delivery documents. This must be kept in a safe place and presented to the technician in the event of servicing or when performing the technical safety check (STK).

Please note that a medical device book is not required in every country.

Note the local requirements and laws.

5.5 Device Parts and Accessories

5.5.1 The device



The Cobra consists of a slit lamp with a built-in Q-Las laser (1) mounted on an instrument table in which a CITO 532 laser is integrated (2) and a display (3). The rotation knob to switch back and forth between CITO and Q-LAS. This is the major function of the device.(4)

5.5.1.1 Instrument table

The instrument table on which the slit lamp is mounted is height adjustable. This allows you to adjust the optimal working height comfortably for each patient. Please note that the table can only take an additional load of 2.0 kg. This means that no objects heavier than 2.0 kg in total may be placed on the table. Also, no persons may sit on the table or support themselves with their entire body weight on it.

If the table is tilted, you can calibrate it. To calibrate the table, move it all the way down and press and hold the button (down arrow) until it comes to a complete stop. After calibration, the table is horizontal again.

5.5.2 Basic equipment

Additionally to the device itself, the basic parts consists of functional parts that are necessary for the operation of the device, as well as additional equipment for safety and functionality.

Necessary parts for operation:

Part	Description	Article number
Mains cable*	Type E/F (EU)	KB13003
	Type G (USA)	KB13004
	Type B cable	KB13005
Interlock plug	Interlock plug for the door interlock	UG01127
Chin rest paper**	Chin rest paper for the patient	SL01115
User manual	User manual	

*The type of the mains cable depends on the destination country. One cable is included in the standard scope of delivery.

**Only chin rest paper provided by A.R.C. laser is allowed to be used with the Cobra. When the enclosed chin rest paper is used up, refills can be obtained from A.R.C. laser.

Additional equipment included:

The packaging equipment consists of:

Part	Description	Article number
Protection carton	Protection carton for the CITO 532 laser box	VP09001
Transport carton	Shipping box incl. distance foam	VP06028
Transport box***	Wooden shipping box	VP01081

***Lasers that are not delivered via freight transportation may not include a wooden shipping box.

other equipment:

Part	Description	Article number
Spare LED Module	Spare module for the LED	BG06565
Dust protection cover	Dust protection cover for the device	ZU01046
Focus rod	Focus rod for adjustment of binoculars	ME06883
Cover for toothed racks	Cover for toothed racks to adjust the slit lamp position	ME06544

5.5.3 Optional accessories

Optional accessories are marketed separately. The products listed can, but do not have to, be purchased from A.R.C. laser.

A.R.C. laser offers **laser safety goggles** for eye safety of additional people in the operating room.

Part	Description	Article-No.
Laser safety glasses (CITO use)	Protection class D LB6+ I LB7 + RM LB8 / Optical density OD8+, 532 nm	AS01003
Laser safety glasses (Q-LAS use)	Protection class D LB 6 + IR LB8 + M LB10 / optical density OD10+, 1064 nm	AS01033

When buying laser safety goggles, please see the requirements in chapter 4.2.

Contact glasses* required for SLT, Iridectomy or capsulotomy procedures can also be obtained from A.R.C. laser:

Application	Description**	Article-No.
SLT	Ocular Latina SLT Gonio Laser Lens	WE01356
Iridectomy	Ocular Abraham Iridectomy YAG Laser Lens	WE01243
Capsulotomy	Ocular Abraham Capsulotomy YAG Laser Lens	WE01244

*All offered contact glass forms are well established for their respective treatment. For further information about the medical treatments or the recommended devices consult the application

manual or your medical device consultant. When buying contact glasses for use with Cobra, make sure that they are suitable for laser treatment.

**The listed designations are the trade names by ocular instruments. They may differ for different providers.

For the personal comfort of the physician and the stabilization of the arm while holding a contact glass, any comfortable **elbow rest** can be used. A.R.C. offers a stackable elbow rest made of foam. (ZU11056, diameter 116mm, height 28mm).

ATTENTION

Only spare parts which are approved by A.R.C. Laser GmbH are to be used with the device. Accessories that have not been approved can significantly impair the safety and reliability of the device.


The use of accessories, transduce and services other than those which the A.R.C. Laser GmbH has determined or provided, may result in increased electromagnetic interference or reduced electromagnetic immunity of the device and lead to incorrect operation.

6 OPERATION

This part of the manual only describes the technical application of the device, without going into the medical application. You will receive separate application manuals from A.R.C. Laser GmbH.

Settings and adjustments should only be made in accordance with the operating instructions. Changes or settings that are not described in this manual can lead to malfunctions.

The laser is in the READY mode during treatment. If you interrupt, the laser must be returned to STANDBY mode. The device must always be switched off when unattended to prevent operation by an unauthorized third party.



The main switch under the table switches the entire system on / off.

ATTENTION

Since the aiming beam takes the same way through the laser transmission system as the working beam, it is a good way to check the integrity of the laser transmission system. If the aiming beam does not appear at the distal end of the laser transmission system, its intensity is weak or if it looks diffuse, this is a possible indication of a damaged or malfunctioning laser transmission system.

ATTENTION

The use of the controls or adjustments for performing procedures other than those specified herein may result in hazardous laser radiation exposure.

6.1 Start the Cobra



Please check the following key points when preparing the device for operation, in order to avoid unnecessary troubleshooting or even malfunctions.

- Make sure the power cord is plugged in.
- Make sure the external door interlock switch is connected.
- Check that the Laser-Stop button is not pressed. If it is pressed, pull it up.
- Turn the rotary knob for device selection to the neutral position (No laser device is selected)
- Activate the main switch at the table
- Switch on the Cobra via the remote-control Key Switch.

Make sure that a sufficient number of protective goggles are ready to be used.



Position Key Switch:



ATTENTION

The Cobra should not be operated for longer than 5 hours at a time. Restart the laser after 5 hours at the latest so that it can carry out the system check during the start routine. This ensures that undetected errors cannot occur.

When the rotary knob on the side of the slit lamp is in the neutral position, i.e., neither the CITO 532 nor the Q-Las is selected, the laser starts normally, and the following image will appear on the display:



The Cobra will now begin to carry out system checks. The system checks for the CITO 532 will start automatically whether the rotary knob is in neutral position or selects the usage of the CITO 532.

To start the initialization process of the Q-Las, turn the rotary knob upwards to select Q-Las.

The display of the slit lamp begins a countdown from 100 to 0. This will take approx. 60 seconds. If the operating temperature has not been reached after this time, the device will start a new warm-up phase and a new countdown will start. After a total of up to three warm up circles, (up to 180 seconds) the device has reached its operating temperature. The LEDs at the slit lamp start to blink - - -

For the Cobra variant with a built-in Q-Las diode version, the following steps are not necessary, and your device is ready to use when the display shows the currently set energy.

If this is not the case:

- **Set the energy level to 5.0 mJ**
- The laser performs 3 simulation shots
- The yellow READY LED flashes
- The Counter will show *0 0 0*
- The current output (5 mJ) is shown in the upper part of the display
- **Now reduce the energy and set your treatment value.**
The laser cannot be switched to READY with 5.0 mJ.

Now Cobra is ready for operation.

6.1.1 Select Application

CITO 532:

To select the CITO 532, turn the rotary knob located on the side of the slit lamp counter clockwise (down) until it clicks into place. The following image will appear on the display:




The main screen of the CITO 532 is then automatically displayed. Use the CITO as described in chapter 6.2 f.

Q-Las:

To select the Q-Las, turn the rotary knob located on the side of the slit lamp in the clockwise direction (up) until it clicks into place. The following image will appear on the display:




Access the system settings for display brightness, sound and language via the gearwheel icon.





Further operation of the Q-Las is described in Chapter 6.3 f.

6.2 Operation of CITO 532



6.2.1 Select doctor / Add doctor

Since it is possible to store multiple users, select a desired user by clicking on the doctor's name at the upper left of the control panel and it will automatically take you to the program stored for that selection. When logging in for the first time, it is necessary to click "New Surgeon" to create a new user.



As soon as several surgeons are entered, they appear in a selection list.



Surgeon List


Select the desired name from the list here.

 EDITING a surgeon

 DELETE from List

A security question prevents accidental deletion.

6.2.2 Selection of Programs




Both programs P01 and P02 can be individually defined and saved.

First set your desired energy and then press the disk symbol to save. The set energy value is now saved. By default, P01 is set to 0.2 mJ and P02 to 0.4 mJ.

Click P01 or P02 selects the respective program.

6.2.3 Energy Configuration

With Plus/Minus the energy can be selected within small steps. Selection using the slider is suitable for larger steps or faster selection.



The following energy values are available:

0.2 till 1.4 mJ in 0.1 mJ steps

1.4 till 2 mJ in 0.2 mJ steps In READY mode, a recalibration is carried out at a value of 1.4 mJ so that the Ready button flashes once.

6.2.4 Laser STANDBY / Ready

This button switches from STANDBY to READY and back. A safety pause of 2 seconds is programmed before the first laser shot. The yellow READY field flashes three times




ATTENTION


If the laser is in the READY mode and is not used for more than 2 minutes, the laser automatically switches to the STANDBY mode.

6.2.5 Counter


The number of shots fired, and the accumulated energy are shown in the lower part of the display (shot counter). To reset the counter to 0, press the sigma symbol (for at least 1s). This will reset the number of shots and the accumulated energy.



You can also access the recorded treatment data via the clipboard symbol (laser log). The system shows the serial number of the device and the average energy set. The total laser shots and the total energy output are displayed here. With the reset button, all counter readings in this and the main mask are reset.



6.2.6 Settings / Submenu



Parameters are set using the up/down arrows



6.2.6.1 Select Language

6.2.6.2 While using the up/down arrows the system language can be selected: DE for German, EN for English Reminder of the Doctor Programs set



and

The system remembers the last selected setting when the paper clip is shown in white. The memory function applies to the doctor and treatment program. At the next start, the system automatically jumps to the last set routine. If the paper clip is dark, nothing is saved or preset.

6.2.7 Release Laser Beam

To trigger the laser, press the trigger on the vertical drive (joystick) of the laser slit lamp. To do this, the laser must be switched to the READY mode.

When a laser shot is fired, a simultaneous warning signal sounds and the red LED on the emergency stop button lights up.



6.2.8 Treatment

Position the patient's chin on the chin rest. Make sure that the patient only touches the application parts. The application parts are the chin rest including handles. Make sure that you and the patient do not touch any other parts at the same time.

Adjust the position of the slit lamp, slit lighting, focus and contact glass as required. Choose your initial values for the laser energy. It is always good to start at a low level and then increase energy to the target during treatment.

Fix the patient's eye with a contact glass and focus the target beam into the eye to be treated.

Select the necessary beam intensity. At this point, you should see that the target steel point is in focus with the microscope.

Press the READY button: the yellow LED flashes for 2 seconds and then changes to steady light. Treatment can start.

The software of the CITO 532 allows the energy delivered to be added up, which is permanently shown in the lower area of the display.

Change of Patient

At the end of a treatment, press the READY button, so you can document the total energy delivered. The trigger does not work in this position.

ATTENTION

Any serious incident that occurs with this laser must be reported to the A.R.C. Laser GmbH and the responsible state authority.

6.3 Operation of the Q-Las



The display shows different information depending on the state of the unit.

6.3.1 Energy Selection


Use the energy selection wheel to change the values. The available energy values can be continuously adjusted from 0.5 mJ to 10 mJ (BurstMode 1), 20 mJ (BurstMode 2) up to a maximum of 30mJ (BurstMode 3).

6.3.2 Burst-Mode (Number of pulses = pulses per shot)

The number of pulses can be varied between 1, 2 or 3. The total energy output (mJ) changes depending on the number of pulses.

By changing the burst mode, the total energy changes as follows:

Energy (mJ) according to burst mode (1,2 or 3 pulses)	
1 Pulse	0.5 mJ to max; in 0.1 mJ steps
2 Pulses	1.0 mJ to max; in 0.2 mJ steps
3 Pulses	1.5 mJ to max; in 0.3 mJ steps



6.3.3 Information LED Energy level (only Q-Las diode version)

This information LED is off as long as the laser energy is exposed within the standard tolerance range of +/- 20 %. As soon as the tolerance of the exposed laser energy is between + 20 % up to + 30 % the information LED lights up in blue and in addition an acoustic signal is emitted to inform the user that the laser is working in the extended tolerance range.



In case that the energy exceeds this extended tolerance range, the laser goes into the test mode. After the test mode has been successfully completed, the treatment can be continued. In case the test mode is not successfully completed a restart of the device is needed.

6.3.4 RESET Counter

To reset the counter to 0, pulse the Reset-button on the remote control: 000 is displayed again.



6.3.5 Laser STANDBY / READY

This button switches from STANDBY to READY and back. There is a safety break before the first laser shot programmed for 2 seconds. During this safety break, the yellow LEDs flashes.

Yellow LED **lights up**: Laser READY

Yellow LED **does not light up**: Laser STANDBY




ATTENTION

If the laser is in READY mode and is not used for more than 3 minutes, the laser automatically switches to STANDBY mode.

6.3.6 Cumulated applied energy

By default the sum of emitted laser shots is displayed. It is also possible to display the applied energy. by pressing the READY LED. The cumulated energy in Joule is displayed for a short time (2-3 s) and the devices changes to STANDBY. The total energy is accumulated until the RESET button is pressed. It resets both the shot counter and the cumulated applied energy.




Note: Your service technician can activate or deactivate this feature.

6.3.7 DualSpot aiming beam

The Q-LAS is equipped with a DualSpot aiming beam as standard. The aiming beam is only visible in READY mode or while setting its brightness in STANDBY mode.

The two beams intersect in the focal plane to form a single spot. The laser exerts its effect on the focal plane. An exception is the focus shift (see chapter 6.3.8).

6.3.8 Adjustments at the laser head



Adjustments at the laser head

6.3.9 Trigger Button at the Joystick

The aiming beam is only visible in READY mode. To trigger the laser, press the trigger on the elevator (joystick). To do this, the laser must be switched to READY mode and the lighting arm must be swiveled out approx. 12 ° to the right or left. An acoustic signal and the display (POS) warns that the arm is not swung out. The laser shot can only be triggered when this is swiveled out far enough. When a laser shot is fired, a warning signal sounds and the LED next to the joystick lights up.



6.3.10 Treatment

Position the patient's chin on the chin rest. Make sure that the patient only touches the application parts. The application parts are the chin rest including handles. Make sure that you and the patient do not touch any other parts at the same time.

Adjust the position of the slit lamp, illumination and focus so that the visible aiming beam points merge into one point.

Choose your initial values for energy and burst mode. It is always good to start with a low energy level and then increase the energy during treatment.

Fix the patient's eye with a contact glass and focus the target beam into the eye to be treated.

Select the necessary beam intensity. At this point, you should see that the target steel point is in focus with the microscope.

Press the READY button: the yellow LED flashes for 2 seconds and then changes to steady light.

Treatment can start.



After setting the first treatment values, it is advisable to first set a laser pulse in the periphery of the retina, or the treatment area, in order to observe the laser effect. Under no circumstances should the first laser pulse be placed centrally or close to the macula.

Change of Patient

At the end of a treatment, press the READY button, so you can document the total energy delivered. The trigger does not work in this position.

ATTENTION

Any serious incident that occurs with this laser must be reported to the A.R.C. Laser GmbH and the responsible state authority.

6.4 Switch Off your Cobra

- 1) Turn the key switch (see chapter 6.1) to the "OFF" position to switch off the laser.
- 2) Remove the key.
- 3) Switch off the device with the main switch on the table.

ATTENTION

DO NOT FORGET TO TURN OFF THE SYSTEM. It should be covered with a dust cover to avoid dust deposits on the optics.

6.5 Laser-Stop

Restart the device after pressing the laser stop (red button):

- 1) Turn the key switch to position „OFF“
- 2) Pull out the red knob on the right side
- 3) Start the device again

The counter is reset automatically. The energy display shows the selected energy value.



ATTENTION

The **Laser-Stop** should only be used in emergencies.

It is located next to the key switch.

7 SPECIFICATIONS

7.1 General

Cooling:	Air internal, active air cooling for CITO 532
Weight:	ca. 51 kg with table and slit lamp
Dimensions:	H 130 cm/B 69.2 cm/T 58.0 cm with slit lamp, height adjustable table 72-94 cm
Control options:	7-segment display, front panel (Q-Las) Digital touch screen
Power transmission:	Coupled with a slit lamp

7.2 Laser data

CITO 532

Laser type	Frequency-doubled Q-switched Nd:YAG Laser
Output	Single Pulse, maximum manual repetition rate: 10 Hz, pulse duration <9 ns

Wavelength	Energy	Tolerance
532 nm	2 mJ max. Settings: 0.1 mJ Steps from 0.2 to 1.4 mJ 0.2 mJ Steps to 2 mJ	±20% of the shown energy according to IEC 60601-2-22

Q-Las

Laser type	Passive Q-switched Nd:YAG Laser
Output	Single Pulse, maximum manual repetition rate: 3 Hz, 2.4 Hz or 1.8 Hz; pulse duration < 5 ns

Wavelength	Energy	Tolerance
1064 nm	40 mJ max. Settings: 0.1 mJ Steps from 0.5 to 10 mJ BurstMode: 1, 2 or 3 Pulses	±20% of the shown energy according to IEC 60601-2-22 extended tolerance range 20 to 30 % of the shown energy according to IEC 60601-2-22 Abs.201.13.2.101 (for Q-Las diode version)

Aiming Beam

Lasertype:	Diode laser
Wavelength:	635 nm, red
Max. output power:	<1 mW
Brightness:	Variable
Mode:	CW, with CITO 532 additional pulsed mode (2 HZ)

7.3 Electrical Connection Data

Power supply connection values 100-240 V, 50Hz/60 Hz, 400 VA

7.4 Classification

Laser class Laser Beam 3B
(Classification EN 60825-1)

Laser class Aiming Beam 2
(Classification EN 60825-1)

Classification according to MDR IIb, Regel 9

Electrical protection class II
(Classification IEC 60601-1)

Certification CE 0123

8 SERVICE

8.1 Introduction

The device was designed, developed and tested according to the latest technical knowledge. We have set the product life to 7 years. In addition, the availability of spare parts is guaranteed by us within a period of 10 years. However, to ensure that everything works properly, we have made it possible for you to carry out a visual check of the status indicators from the outside.

ATTENTION

There is no need for the laser user to perform routine or service work within the laser system. All adjustments and calibrations that require the protective housing must be carried out by trained service personnel. This also includes cleaning and cleaning the optics within the laser.

8.2 Safety Check (STK)

Once within 24 months, the laser must be subjected to a safety check (STK) by trained personnel. The execution of the STK and any faults are to be noted in the device book.

Please note that a medical device book and a regular technical safety check are not required in every country. Note the local requirements and laws.

Scope of Safety Check

Visual Check

- Laser marking (laser class, max. Power, wavelength)
- Information signs/warning signs; properly and completely attached
- Instructions for use/medical device book
- Condition of the supply line
- Condition of the goggles/protective device
- Overall condition

Functionality Check

- Footswitch/joystick
- Beam guidance system/coupling/decoupling/pilot laser
- Check operating foil (touch screen)

Check the necessary monitoring/safety display and signaling unit

- Laser safety glasses
- Emission controls (acoustic, visual)
- Power meter (comparison internal, external)
- Key switch
- Laser stop switch (check for function)
- Interlock device (check for function)

Electrical Safety Check

- According to IEC 62353 or DIN VDE 0751, Part 1

Output power measurements

- Check the set power with an external power meter (permissible tolerance + 20%)

ATTENTION

If one or more safety-related points are objected to after the safety-related inspection (STK), the device may no longer be operated.

8.3 Care by the user

The following maintenance instructions can be carried out by the user. These serve to make your work easier. The system must be disconnected from the mains for cleaning. Use a damp, but never wet, soft cloth.

It is recommended to disinfect all parts touched by patients in order to avoid the transmission of pathogens. Handles, headgear and chin rest can be disinfected at any time. Chinrest paper is used to further improve hygiene. The top sheet can be easily removed after use or before each new treatment. DescoseptAF * or a comparable disinfectant (e.g. Mikrozid AF) can be used for surface wiping disinfection. (* DescoseptAF: Dr. Schumacher GmbH (www.schumacher-online.de)) DescoseptAF * solution contains approx. 42% ethanol and approx. 0.05% didecyldimethylammonium chloride.

Other disinfectants can be used, as long as they are not aggressive or contain acids. For example, Agents based on quaternary ammonium compounds such as TPH protect (from Schülke) or Mikrobac® forte from Bode.

The manufacturer information for combined devices - for example slit lamps - must be observed!

ATTENTION

During cleaning, the device must always be switched off. Water or disinfectant should never be applied directly to the surface of the device. Make sure to dabble the cloth to then clean or disinfect the device parts.

8.3.1 Inspection and cleaning of external optics

Check the accessible surfaces of the optics for possible contamination. The surfaces of the slit lamp can be cleaned with a soft cloth and distilled water or a mixture of distilled water and approx. 10 % alcohol.

8.3.2 Inspection and cleaning of internal slit lamp optics

Regularly check all optics on the slit lamp - especially the plug-in mirror and the lens above it - for contamination.

A contaminated or dirty lens could lead to a reduction of energy when the laser beam passes through. The following problems can result:

- The set performance is no longer achieved during the test and calibration.
- A change in the beam behavior and the associated active beam during the operation.

Remedial Action:

- 1) Use non-alcoholic, not soaked cotton swabs, and try to clean the outside of the lens with distilled water. If this does not lead to success, please continue with step number 2.
- 2) Take a few drops of methanol and gently drizzle the end of the cotton swab, then do as in step number 1.

ATTENTION

Always wipe the optical surface in one direction. Never go back in the other direction with the same cotton swab. Any particles that have already settled in the cotton swab would scratch the optical surface again when moving back and thus cause irreparable damage to the optics.

8.4 Slit lamp


You can find detailed information on the care and maintenance of the slit lamp in the accompanying manual. It is essential to follow these instructions.

8.4.1 Dust Cover

After each use, the dust cover and the red display protective cover should be replaced to keep the surfaces clean.

8.4.2 Testing the Slit Lamp Focus

1. Use the focus rod and insert it into the hole in the slit lamp axis. Follow the slit lamp manual that comes with the device.
2. Turn on the laser. Make sure that the laser is not triggered during this time. The laser must be in the READY mode. Center the aiming beam in the middle of the visual field on the focus rod.
3. Check that the aiming beam is circular and symmetrical.
4. A disturbed target beam can be an indication of a faulty system; thus the aiming beam can be a tool to check the integrity of the system.







ATTENTION

Press the eyepieces gently into the holder so that they are held completely in the binocular piece. Adjust each eyepiece to your satisfaction so that you can see a sharp image of the focus rod through each eyepiece.

8.5 Troubleshooting

8.5.1 In case of an error, warnings will be displayed in the Touch Screen (CITO 532)

 <p>Error! Code number: F06</p> <p>An error has occurred. Please restart the laser. If the error persists, please call your service provider.</p> <p>A.R.C. Laser GmbH Phone: +49-911-217790</p>	<p>F-errors are hardware errors Common F-errors:</p>	
	<p>F01</p>	<p>F06, F07, F08</p>
 <p>Interlock message!</p> <p>Door of the operation room is open or the interlock plug is not connected.</p>	<p>Sub-performance Please restart the system.</p>	<p>Problems with the trigger. Make sure that the trigger is not pressed continuously. Please restart the system if necessary.</p>
 <p>Interlock message!</p> <p>Foot switch / Manual trigger pressed! Please release the switch before entering READY-mode.</p>	<p>If the error cannot be remedied, contact the service of A.R.C. Laser GmbH.</p>	
 <p>Interlock message!</p> <p>Trigger was not completely activated.</p>	<p>Door contact error Close the door to the treatment room or check the cable connection to the door contact switch.</p>	
	<p>Joystick or footswitch pressed. Release the switch to be able to continue working.</p>	
	<p>Joystick or footswitch only half-pressed. Release the switch to be able to continue working</p>	

<p>Interlock message!</p>  <p>Ambient temperature too low! Risk of condensation.</p> 	<p>The ambient temperature is too low, the laser should be warmed up to continue working. There is a risk of damage to the laser.</p>
<p>Interlock message!</p>  <p>Ambient temperature too high. No cooling of the device possible.</p> 	<p>The ambient temperature is too high, the laser should be cooled in order to continue working.</p>
<p>Interlock message!</p>  <p>Internal test of the output power in progress. Please wait...</p>	<p>Please wait</p>
 <p>The service schedule for your laser has exceeded! Please contact your service provider: A.R.C. Laser GmbH Phone: +49-911-217790</p> <p>Continue</p>	<p>Contact Service of A.R.C. Laser GmbH</p>

8.5.2 In case of an error, warnings will be displayed in the 7-segment display (Q-Las)

Possible Error		Messages and solutions	
F00	Watchdog error	Device does not proceed past system check	Hardware Error Please restart the system. If the error occurs regular, contact your local sales & service-representative of A.R.C. Laser
F01	Reduced output -20%	Permanently checked by test shots during the system start-up and during the operation	Hardware error Please restart the system
F02	Higher output +20% (+30% Diode version)	Permanently checked by test shots during the system start-up and during the operation	Hardware error Please restart the system
F04	Safety Shutter error (Timeout or open/close failure)	Permanently checked by test shots during the system start-up and during the ready mode	Please contact your local sales & service-representative of A.R.C. Laser
F06 F07 F08	Foot switch short circuit – pre / post	Checked during the system start	Exchange the footswitch. If problem persists, please contact your local sales & service-representative of A.R.C. Laser
F11	Unable to set power - pol cube error	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F21	Error READY button (short circuit)	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F38	Error DA converter	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F44 Flash lamp version	Charging regulator of power supply without function	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F44 Diode version	Cavity temp. or diode driver/pulse monitoring error	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F47	Performance check, Checksum error	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F50	I2C-bus error / memory error	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F51	UI-button defective	User interface error at the front panel	Please contact your local sales & service-representative of A.R.C. Laser
F58	Feedback voltage to high (Flash lamp version) or diode current out of range (Diode version))	Checking the maximum voltage at the internal charge regulator (Flash lamp version), Checking if diode current is to high or to low (Diode version),	Please contact your local sales & service-representative of A.R.C. Laser
F61	Low voltages at the charging regulator	Checking the minimum voltage at the charging regulator	Please contact your local sales & service-representative of A.R.C. Laser
F62 Flash lamp version	Energy test error or burst mode voltage monitoring error	Checking if the energy or the burst voltage monitoring is function properly (d)	Please contact your local sales & service-representative of A.R.C. Laser

F62 Diode version	Light measurement unit error	Test of light measurement unit at system check	Please contact your local sales & service-representative of A.R.C. Laser
F63	Temperature/humidity error	Checking the temperature, humidity respectively	Decrease the temperature, humidity respectively.
F65 Diode version only	Head-up display connection lost	Checking the connection to head up display during system check	Please contact your local sales & service-representative of A.R.C. Laser
FOC	Position of the focus shift dial	Position of the focus shift is constantly checked	Check the position of the focus shift. It should be locked.
IL1	Internal Temperature	The internal temperature is out of permitted range. The laser should be allowed to cool down.	Please press the Reset-Button to delete the message.
IL2	Time Delay Trigger	Pressing the trigger only half way down or too slowly may cause this error message.	Please press the Reset-Button to delete the message. You immediately can proceed working afterwards.
IL3	Internal power monitoring	Error when changing the energy while the device changes from STANDBY to READY	Device goes to STANDBY. Press the counter reset button to clear the display.
A or POS	Position of the arm: PCL5-SH and Z-version	The position of the lighting arm is constantly checked	Swing out the arm (or for Diode version: lower prism out of beam path)
Possible Error		Messages and solutions	
F00	Watchdog error	Device is hanging during system check	Hardware Error Please restart the system. If the error occurs regular, contact your local sales & service-representative of A.R.C. Laser
F01	Reduced output -20%	Permanently checked by test shots during the system start-up and during the operation	Hardware error Please restart the system
F02	Higher output +20%	Permanently checked by test shots during the system start-up and during the operation	Hardware error Please restart the system
F04	Optokoppler Safety Shutter Time Out	Permanently checked by test shots during the system start-up and during the ready phase	Please contact your local sales & service-representative of A.R.C. Laser
F06 F07 F08	Foot switch short circuit – pre / post	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F11	Unable to set power	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F21	Error READY button	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F38	DA converter	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F44	Charging regulator of power supply without function	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F47	Performance check, Checksum error	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser
F50	I2C-error	Checked during the system start	Please contact your local sales & service-representative of A.R.C. Laser

F51	Keys error	User interface error at the front panel	Please contact your local sales & service-representative of A.R.C. Laser
F58	Pulse Threshold Voltage	Checking the maximum charge voltage at the internal charge controller	Please contact your local sales & service-representative of A.R.C. Laser
F61	Low voltages at the charging regulator	Checking the permissible voltage at the charging regulator	Please contact your local sales & service-representative of A.R.C. Laser
F62	Energy test error or burst mode voltage monitoring error	Checking if the energy or the burst voltage monitoring is function properly (d)	Please contact your local sales & service-representative of A.R.C. Laser
F63	Temperature/humidity error	Checking the temperature, humidity respectively	Decrease the temperature, humidity respectively.
F0C	Position of the focus shift dial	Position of the focus shift is constantly checked	Check the position of the focus shift. It should be locked.
IL1	Internal Temperature	The internal temperature is out of permitted range. The laser should be allowed to cool down.	Please press the Reset-Button to delete the message.
IL2	Time Delay Trigger	Pressing the trigger only half way down or too slowly may cause this error message.	Please press the Reset-Button to delete the message. You immediately can proceed working afterwards.
IL3	Internal power monitoring	Error when changing the energy while the device changes from STANDBY to READY	Device goes to STANDBY. Press the counter reset button to clear the display.
A or POS	Position of the arm: PCL5-SH and Z-version	The position of the lighting arm is constantly checked	Swing out the arm

8.6 Disposal

The relevant, locally applicable laws and regulations must be observed during disposal. Under no circumstances should the device be disposed of with domestic waste.

A.R.C. Laser GmbH is happy to assist with the disposal of the system; Costs for the proper return of the laser to A.R.C. Laser GmbH is the responsibility of the buyer. Please contact our service department.



9 CUSTOMER SERVICE

9.1 Warranty Information

A.R.C. Laser GmbH grants you a two-year guarantee. Parts that have a defect will be replaced free of charge within two years. All add-on and purchased parts are exempt from this guarantee. Our guarantee extends to the repair or replacement of defective parts. However, we reserve the right to renew entire assemblies and adapt them to technical progress.

Repairs by third parties or changes to the device will void the warranty. The use of other parts that have not been accepted with the device or obtained from other suppliers will also void the warranty. The attachment of parts or assemblies or other changes to the device requires the express written confirmation by A.R.C. Laser GmbH.

9.2 Warranty, Shipment, Packing

A warranty claim for defective parts, malfunction or damage to the housing of the device must be submitted to A.R.C. Laser GmbH within 24 hours. Parts that are returned during the warranty period (at the express request of A.R.C. Laser GmbH) must be confirmed in writing by A.R.C. Laser GmbH. Detailed packaging instructions and information on how to return the device are provided by A.R.C. Laser GmbH. The return must be insured and paid for by shipper. The insurance and transportation costs are not covered by A.R.C. Laser GmbH. The choice of the return is made by the A.R.C. Laser GmbH communicated to the customer. Changes and amendments in the carrier or the shipping method can lead to delays in transport and processing. All components to be changed under the warranty claim are manufactured by A.R.C. Laser GmbH renewed free of charge within the guarantee period. We reserve the right to make changes to the design of the device - if it appears necessary - to increase the safety or the functionality of the device. The responsibility for the design as well as for changes in the device lies solely with A.R.C. Laser GmbH. Changes will be communicated to the customer and accordingly carried out at A.R.C. Laser GmbH.

9.3 Sales and Service Information

For sales and service information, please contact A.R.C. Laser GmbH or our local distributor.

10 GUIDELINES AND MANUFACTURERS DECLARATION

10.1 Electromagnetic Emissions		
The laser is intended for use in an environment as specified below. The customer or user of the laser should ensure that it is operated in such an environment.		
Immunity tests	Compliance	Electromagnetic environment - guideline
RF-Emissions CISPR 11 150 kHz – 30 MHz	EN 55011 Group 1/Class B	The laser is only suitable for the environment in professional healthcare facilities.
RF- Emissions CISPR 11	EN 55011 Group 1/Class B	
Harmonic emission	IEC 61000-3-2 Class A	
Voltage fluctuations/flicker	IEC 61000-3-3	

10.2 Electromagnetic Immunity (1)


The laser is intended for use in the electromagnetic environment specified below. The customer or the user of the laser should ensure that it is used in such an environment.

Immunity tests	IEC 60601-Test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic discharge (ESD)	± 2 kV, ± 4 kV, ± 6, ± 8 contact discharge; ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air discharge	IEC 61000-4-2	Floors should be made of wood, ceramic or stone. If the floor is covered with a synthetic material, the relative air humidity should be at least 30%.
Electrical fast transient /burst	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz	IEC 61000-4-3	Only the voltage of a typical professional healthcare facility may be used.
Radiofrequency electromagnetic fields in the immediate vicinity of wireless communication devices	3 V/m 80 MHz to 2,7 GHz (see Table 10.4)	IEC 61000-4-3	Only the voltage of a typical professional healthcare facility may be used.
Magnetic fields with energetic design frequencies	30 A/m 50 Hz or 60 Hz	IEC 61000-4-8	Magnetic fields at the grid frequency that are usually available in professional healthcare facilities may be used.
Electrical fast transient /burst	± 2 kV for power lines ± 1 kV for IO-lines 100 kHz repetition frequency	IEC 61000-4-4	Only the voltage of a typical professional healthcare facility may be used.
Surge voltages (Surges), Line against line	± 0.5 kV, ± 1 kV	IEC 61000-4-5	Only the voltage of a typical professional healthcare facility may be used.
Surge voltages (Surges), Line against grounding	± 0.5 kV, ± 1 kV, ± 2 kV	IEC 61000-4-5	Only the voltage of a typical professional healthcare facility may be used.
Conducted disturbance variables, induced by high-frequency fields	3 V 0.15 MHz to 80 MHz 6 V in ISM-frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	IEC 61000-4-6	Only the voltage of a typical professional healthcare facility may be used.
Voltage dips	0 % UT; ½ cycle at 0, 45, 90, 135, 180, 225, 270 and 315 level 0 % UT; 1 cycle at 0 and 180 level And 70 % UT; 25/30 cycles at 0 and 180 level	IEC 61000-4-11	Only the voltage of a typical professional healthcare facility may be used. It is recommended to use an uninterruptible power supply.
Power interruption	0% UT; 250/300 cycles at 0 and 180 level	IEC 61000-4-11	Only the voltage of a typical professional healthcare facility may be used. It is recommended to use an uninterruptible power supply.

NOTE: UT is the AC mains voltage prior to application of the test level.

10.3 Electromagnetic Immunity (2)

The laser is intended for use in the electromagnetic environment specified below. The customer or the user of the laser should ensure that it is used in such an environment.

Immunity tests	IEC 60601-Test level	EMV standard	Electromagnetic environment – guidelines
Conducted disturbance variables, induced by high-frequency fields	3 V 150 kHz to 80 MHz	IEC 61000-4-6	<p>Portable and mobile RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the CLASSIC 514, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.</p> <p>The field strength of stationary radio transmitters is, as determined by an electromagnetic site survey, at all frequencies smaller than the compliance level.</p>
Radiofrequency electromagnetic fields in the immediate vicinity of wireless communication devices	3 V/m 80 MHz to 2.7 GHz	IEC 61000-4-3	<p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div> <p>If higher IMMUNITY TEST LEVELS than those specified in Table 9 (IEC 60601-1-2) are used, the minimum separation distance may be lowered. Lower minimum separation distances shall be calculated using the equation specified in 8.10 (IEC 60601-1-2).</p>

10.4 Recommended separation distances between portable and mobile RF telecommunications equipment and the laser

The laser is intended for use in an electromagnetic environment in which the RF disturbances are controlled. The customer or the user of the laser can help to avoid electromagnetic interference by maintaining the minimum distance between portable and mobile HF telecommunication devices (transmitters) and the laser - depending on the output power of the communication device, as stated below.

Test frequency	Frequency band	Radio service	Modulation	Maximum performance	Distance	Immunity test level
MHz	MHz			W	m	V/m
385	380 to 390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
450	430 to 470	GMRS 460 FRS 460	FM ± 5 kHz Hub 1 kHz Sinus	2	0.3	28
710 745 780	704 to 787	LTE Band 13,17	Pulse modulation 217 Hz	0.2	0.3	9
810 870 930	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
1720 1845 1970	1700 to 1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1,3,4,25, UMTS	Pulse modulation 217 Hz	2	0.3	28
2450	2400 to 2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240 5500 5785	5100 to 5800	WLAN 801.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9

Notes:

A.R.C.

L A S E R

made in Germany ● ● ●

A.R.C. Laser GmbH

Bessemerstr. 14

D-90411 Nürnberg

Tel: +49 911 21779-0

Fax: +49 911 21779-99

www.arclaser.de