

Q-Las: The Nd:YAG in your daily routine for Capsulotomy and Iridotomy

Q-LAS the standard laser in your daily routine

A.R.C.
LASER
made in Germany



Ideal for capsulotomy and iridotomy

Q-LAS the standard laser in your daily routine

Advantages

- Perfect optic for anterior segment
- Choose from parallel or convergent binocular
- Dual spot aiming beam
- Clear filter protection

Optional

Illumination tower (prism) lower to keep
0° straight position with the illumination light.



Q-Las YAG-Laser

32mJ max. at 3
pulses



Energy: 0.5mJ to
30 mJ

Focal shift:
150/300µm
posterior

Q-LAS the standard laser in your daily routine

Secondary Cataract
Treatment of patients
suffering from posterior
capsule opacification


opening in the capsule
to remove fibrosis after
cataract surgery to
eliminate opacification



Q-LAS the standard laser in your daily routine

Iridotomy:
Treatment of patients
suffering from angle
closure glaucoma

preventive intervention
for planar anterior
chambers/ narrow
chamber angle




Q-LAS the standard laser in your daily routine

CONVENIENT CONTROL SETTINGS

(Clever details for more comfort)

- Continuous energy levels from 0.5 to 10 mJ
- Prominent buttons, clear, bright display




- Choose your favorite eye piece:
- Either convergent or parallel optics.

1064 nm

Ideal Position of the controller: directly under the binoculars, adjustable angulation


Q-LAS the standard laser in your daily routine




Q-LAS the standard laser in your daily routine

State-of-the-art-capsulotomy laser:

- Perfect optic for anterior segment
- Choose from parallel or convergent binocular
- Dual spot aiming beam
- Head up display (optional)



Q-LAS the standard laser in your daily routine



Power Selection

Use the energy knob to increase or decrease the energy emitted during one shot of the laser.

Focus Shift

The laser focus will be moved posterior:

0: Area of laser effect is $\sim 30 \mu\text{m}$ behind the focal plane

150: Area of laser effect is $150 \mu\text{m}$ behind the focal plane

300: Area of laser effect is $300 \mu\text{m}$ behind the focal plane

A.R.C. Laser only offer the posterior focus shift, as there is no application requires an anterior shift

Aiming Beam

Turn the knob to increase or decrease the brightness of the aiming beam.

Rotation knobs on both sides for left and right handed use

Q-LAS the standard laser in your daily routine

Unique features:

1. Selection of eye piece
 - parallel
 - convergent

parallel



convergent



Q-LAS the standard laser in your daily routine

Unique features:

2. Illumination tower (prism) lower to keep 0° straight position with the illumination light.



Laser can be exposed with illumination in swung in position




Illumination tower need to be in swing out position to laser

Q-LAS the standard laser in your daily routine

Unique features:

2. Head-Up display displays energy information in the microscopic view.





Q-LAS the standard laser in your daily routine

Good to know:

1. Refraction of eye piece is not 0; mostly between -1 and -2 due to the building structure of the laser including filters and head up display
→ Please adjust refraction before use
2. Standard microscope binocular is parallel, it is possible to change to convergent during the order process in case of doctors request
→ Parallel has good focus characteristic, therefor ideal for anterior treatments
→ Convergent has good focal depth, therefore it is ideal for porsterior treatments
3. 2 spot aiming beam instead of 3 spot to reduce the numbers of reflection for a better view and positioning

- Coming soon: Q-Las D – the diode pumped Nd:YAG

Energy stability: Diode vs. flash lamp pumping



Q-LAS-D the latest advance in Nd:YAG laser

Pulse repetition rate:

Flash lamp laser have restrictions in time because of capacitor recharging!!! Most systems: 3 Hz

Diode pumped systems do not need conventional capacitor recharging:

Fast firing possible: (e.g.) 5 Hz (actual techn. spec. will be changing to higher rep. rates)

Q-LAS-D the latest advance in Nd:YAG laser

Further technological advantages:

Reliability, live time:

No high voltage in the laser head! No functional/ electronical problems with humidity.