

COBRA

SLT and YAG combination with 2 separate lasers

COBRA – SLT and YAG combination



SLT and YAG combination

COBRA


- Quick and easy change between both lasers
- safe
- 2 complete separated lasers devices – no overlapping technical parts

COBRA – SLT and YAG combination



- Quick and easy change between both lasers
- Rotation knob to switch from YAG to SLT

COBRA – SLT and YAG combination



- Anterior segment optics (parallel eye piece)
 - convergent optics optional
- 5 times mag. changer: 6.3, 10, 16, 25, 40

COBRA – SLT and YAG combination



Table:

- Wheel chair accessible
- Electronic height adjustment
- Table wheels optional

COBRA Competitor Device

Smart SELECTA Duet Lumenis



Competitor Lumenis:

- One cavity only (one laser system!)
- Parallel optics!?
- 3 Hz (SLT) trigger limited
- 6, 9, 16, 28, 42 mag. changer optics

COBRA Competitor Device

tango
Ellex



Competitor Ellex:

- One cavity only (one laser system!)
- Parallel optics!?
- 3 Hz trigger limited

COBRA Competitor Device

OptoYag&SLT M
OPTOTEK Medical



Competitor Optotek:

- Two cavities
- Parallel optics!?
- 3 Hz (SLT) trigger limited
- 10, 16, 25 x Magnification
- 6, 10, 16, 25, 40 mag. changer optics (optional)



A.R.C. Laser GmbH
Bessemerstraße 14

90411 Nürnberg
Deutschland

