

Device	FOX IV 810	FOX 810
Company	A.R.C. Laser GmbH	A.R.C. Laser GmbH
Picture		
Wavelength	810 nm diode laser	
Indications	CPC μCPC Endocoagulation Retinocoagulation with LIO DCR	
Power Supply	Battery powered, mains supply	
Dimensions	H: 17.4 cm, W: 14.2 cm, D: 16.3 cm	
Power	Max. 8 Watt, safety feature: ant./ post. Treatment up to 3 W	Max. 8 Watt
Voice Output	Yes (10 s 10 J interval) (selection of 5 languages)	No
Counter	Energy and time. Count-down and Count-up mode	Energy Count-up mode
Metronome function	Yes Audio signal every 1 s or 1 J	No
User Interface	Colour touch display (H: 11 mm, W: 6.5 mm) and side bar	Black-and-white touch display (H: 5.5 mm, W: 7.5 mm) and buttons
Colour Codes Applications	Yes	No
Housing	metal	plastic
Weight	1.53 kg	1.2 kg
Power and laser key	Power on/off switch + Password	Power on: key combination + Password Power off by closure of stand
Applicators	Glaucoma Probes, Endo Probes, DCR Probe, LIO	