



**WOLF**  
SURGICAL LASER

**A.R.C.**  
LASER

**PHOTOANGIOLYTIC  
NON-CONTACT CUTTING  
FIBER GUIDED LASER**

Endoscopic  
Office-Based  
Operating Theater

Join the  
**TruBlue** Revolution  
MADE IN GERMANY



# Join the **TruBlue** Revolution

MADE IN GERMANY



## TRULY FLEXIBLE



## TRULY VERSATILE



## TRULY MOBILE



## TRULY NOISELESS



**WOLF**  
SURGICAL LASER

The smallest "CO<sub>2</sub> + KTP"-like laser unit for coagulation, vaporisation and bloodless surgery inside an extremely flexible glass fiber.

## VERSATILITY MEETS EFFICIENCY

The Wolf TruBlue gives you a vast versatility in applications. You are also free to apply the laser in contact or non-contact mode to meet the most beneficial requirements of your desired application in:

**Rhinology, Laryngology,  
General surgery, Otology**


## TRULY CONVINCING

The WOLF TruBlue cuts similar to CO<sub>2</sub> laser in non contact and simultaneously shows superb coagulation characteristics. Thanks to those properties and its high absorption in hemoglobin while being transparent for water (photoangolytic effect) the WOLF TruBlue offers a spectrum of application that seeks its equal.

## TrueBlue

DON'T CALL IT CO<sub>2</sub>

The shockproof and portable WOLF offers unprecedented design, size and function. The laser is supporting a variety of fibers with 300, 400 or 600 micron, providing you with high flexibility. Advances in fiber optics have made them ideally suited for endoscopic use, to travel through routes in the human body where no hand or scalpel has gone before.



**A.R.C.  
LASER**

**A.R.C. Laser GmbH**  
Bessemersstraße 14  
90411 Nuremberg  
Germany

T +49 911 217790  
F +49 911 2177999  
info@arclaser.de  
www.arclaser.de

